

The Right Sonicare Brush Head for Each Sonicare Toothbrush Handle

Each Sonicare toothbrush is designed to work with certain Sonicare brush heads to deliver exceptional results. By choosing the right brush head and changing it every three months, you can ensure that your Sonicare system is performing at its best.



All-Around Clean

ProResults Standard Brush Head

The larger surface area and contoured bristles work with your choice of brushing modes on the toothbrush handle to thoroughly clean teeth and invigorate gums.

Use with: FlexCare+, FlexCare, HealthyWhite



Focused Clean

ProResults Compact Brush Head

Helps you focus the dynamic cleaning action generated by the Sonicare handle on hard-to-reach areas where plaque can accumulate.



All-Around Clean

e-Series Standard Brush Head

e-Series brush heads are angled for better reach within the mouth, providing a thorough, everyday clean while massaging and stimulating gums.

Use with: Xtreme, Essence, Elite



Focused Clean

e-Series Compact Brush Head

The e-series compact brush head helps you focus on hard-to-reach areas where plaque and stains can hide.



Kid-Friendly Clean

Sonicare For Kids Brush Head

Two brush head sizes complement the two power levels on the toothbrush handle to thoroughly and gently clean young teeth and ensure that the brush grows with your children.

Use with: Sonicare For Kids



Change Your Brush Head When the Seasons Change

Bristles exhibit fatigue and brush heads are less effective after three months of normal use, even though these changes are only visible under a microscope. So to keep your Sonicare toothbrush at its best, replace your brush head four times a year.

